

Moving to the FreeStyle Libre 2 system is easy.

Follow the guidance below.

**YOUR PRESCRIPTION HAS BEEN CHANGED FROM
FREESTYLE LIBRE 1 SENSORS to FREESTYLE LIBRE 2 SENSORS.
Follow the steps below to begin using the new Freestyle Libre 2 system**

**CHOOSE WHICH DEVICE YOU WOULD LIKE TO RECEIVE
ALARMS ON?¹**

PHONE



YOU'RE READY TO GO!

Using the FreeStyle LibreLink app, you can start scanning FreeStyle Libre 2 sensors straight away²
Please check that their FreeStyle LibreLink app version is 2.5 or above. See compatibility guide for details.



READER



REQUEST A READER

Go to the Reader Replacement Portal and complete the online request form
FreeStyleLibre.co.uk/replacement

REQUEST VALIDATION

- Enter your current FreeStyle Libre reader serial number
- No requirement to send the current reader back

DELIVERED DIRECTLY TO YOU

A replacement reader is delivered directly to you as the FreeStyle Libre user

Images are for illustrative purposes only. Not actual patient data.

1. Patients choose which device they want to receive alarms: FreeStyle Libre 2 reader or FreeStyle LibreLink app. They must start their FreeStyle Libre 2 sensor with that selected device. Once the patient scans their FreeStyle Libre 2 sensor with that device, they can receive alarms only on that device. 2. The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink requires registration with LibreView.

NOTE: If you are currently using your Freestyle Libre 1 reader to scan your Freestyle Libre sensor, you need to order your free upgraded reader from Abbott AS SOON AS POSSIBLE.

Order on line at: www.freestylelibre.co.uk/replacement OR call customer services on 0800 170 1177. **IMPORTANT:**

- You can only use the Freestyle Libre 2 reader to scan Freestyle Libre 2 sensors.
- You can only use the Freestyle Libre 1 reader to scan Freestyle Libre 1 sensors.